

ings = TALL bucket. so,... no





Select bucket to deploy from new, clean supplies. Note, washers may be needed to raise the bucket when it is deployed.



Trigger sensor by passing fingers through the opening 5 times. Turn off power to collector when it is open.



NTN: N-CON

Sample Change-out (refer to SOP for complete instructions: http://nadp.slh.wisc.edu/lib/manualsSOPs.aspx

Use lid bag to place new, clean lid on the deployed bucket.



Seal lid on bucket using lid bag as a glove.



Remove sealed bucket from collector and set aside.



Deploy new, clean bucket using bucket bag as a glove.

Turn power on to the collector, allowing it to close.

Complete blocks 1-4, and 10 of the field form.

Take the sealed bucket to the field laboratory to weigh and decant the sample for shipment to the CAL for analysis.

NADP OFFICE USE ONLY 1. SITE FIELD BUCK

Revised: 22 June 2018